



Editorial Contact:

Jana Stender
CPC

Tel: +1 952-380-6263

Email: jana.stender@cpcworldwide.com



Hi-res photos available [here](#).
Informative video available [here](#).

New AseptiQuik® Connector For Higher Flow in Bioprocessing

CPC Genderless Connector for High Volume, Large Flow

St. Paul, Minn. – April 17, 2018 – [CPC \(Colder Products Company\)](#), maker of genderless AseptiQuik® sterile connectors, introduces a one-inch option to deliver higher flow rates for high-volume, large-flow transfer lines in full-scale bioprocessing operations. The product introduction coincides with the INTERPHEX conference in New York.

“Today’s bioprocessing market needs higher flow rate connectors,” said John Boehm, CPC Business Director, Bioprocessing. “Customers want simple to use, sterile connectors that ensure flow path sterility and reliability at higher volumes. The new one-inch AseptiQuik L connector is the easiest to use, large-format, genderless connector on the market.”

A one-inch connector enables 3.5 times more flow than a 1/2-inch connector—critical performance for large-volume bioreactors. Genderless connectors simplify use, operator training and inventory management and avoid the need for adapters, which can leak or pose a contamination risk. A robust design provides reliable performance without the need for clamps, fixtures or tube welders.

The AseptiQuik L Series’ integrated pull tabs act as protective covers, reducing complexity and ensuring simultaneous removal of both membranes. The connectors’ easy “Flip-Click-Pull” assembly process and audible “click” confirming connection also help reduce operator error.

For current AseptiQuik users, validation, specification and implementation of the AseptiQuik L connector is seamless. The new connector uses the connection process as other AseptiQuik products. This simplifies training and standard operating procedure (SOP) development. Additionally, a data code on the connector offers easy traceability and immediate access to important product information with a quick scan. The data code connects users to product specifications, date of manufacture, validation reports and training videos.

For more information about the performance and versatility of the genderless AseptiQuik L Series or any of the other 10,000+ innovative connection solutions CPC offers, visit cpcworldwide.com.

About AseptiQuik® Connectors

AseptiQuik® Connectors provide quick and easy sterile connections, even in non-sterile environments—a critical capability for biopharmaceutical and cell and gene therapy manufacturers. Featuring a wide range of options including 1/8- to 1-inch sizes and genderless and gendered connections, the AseptiQuik® connection technology allows easy media transfer with less error risk. Robust, reliable performance eliminates the need for clamps, fixtures or tube welders, delivering sterile, high-quality single-use connections.

###

About CPC

CPC (Colder Products Company) is the leading provider of quick disconnect couplings, fittings and connectors for the life sciences, bioprocessing, industrial and chemical handling markets. CPC is an operating company within Dover Corporation. For a free catalog or more information, contact: Colder Products Company, 1001 Westgate Drive, St. Paul, MN 55114. Phone: 651-645-0091. Toll-free: 800-444-2474. Fax: 651-645-5404. Web: cpcworldwide.com

CPC, Colder Products, Colder Products Company and Dover are registered trademarks with the U.S. Patent & Trademark office.